

We care for every drop



Accessories



IV-Stands

for single pumps or docking stations
height adjustment up to 2.4 m
five antistatic castors and 2 brakes

- IV-stand with one pole
- ARGUS 100 with double pole
- ARGUS 160 with double pole for 2 docking stations

REF. 90.001 IV-Stand

REF. 90.019 ARGUS 100

REF. 90.024 ARGUS 160 (see photo)



ARGUS Docking Stations

compatible with ARGUS infusion and syringe pumps

available in 2 sizes:
with 3 plug-in positions (reference 60)
with 5 plug-in positions (reference 100)

Version P

- power supply, 230 V

Version M

- power supply, 230 V
- RS-232 interface
- interface for Barcode Reader and for staff alert system



REF. 90.202/60 P

REF. 90.200/100 P

115 V on request

REF. 90.102/60 M

REF. 90.100/100 M

115 V on request



Power distributor

- with 4 cables, extendable to a length of 1.80 m
- standard plugs for medical devices
- with universal clamp for poles up to dia. 38 mm

REF. 90.009



Barcode reader with holder

compatible with ARGUS Docking Stations M

- for rail fixation up to dia. 38 mm
- barcode software for safe pump programming on request

REF. 90.151

Barcode reader

- without holder



REF. 90.150

Holder for barcode reader

- spare part to REF. 90.151



REF. 10.130

Drop detector

compatible with ARGUS infusion pumps

- spare part
- integral part of every infusion pump



REF. 10.089



Bottle holder

to be used with ARGUS infusion pumps

- easy insertion
- available in two heights: 45 and 65 cm

REF. 11.005 / 45 cm

REF. 11.043 / 60 cm



Carrying handle

compatible with ARGUS syringe pumps

- anytime mountable
- pump remains docking compatible

REF. 90.050



Table support

to be used with ARGUS syringe pump

- anytime mountable

REF. 90.051



Transport unit

compatible with ARGUS infusion and syringe pumps

- without mains connection
- with 3 plug-ins for 3 syringe pumps or 1 infusion and 1 syringe pump

REF. 90.052



Combi Clamp Standard and for rail fixation

for ARGUS infusions and syringe pumps

REF. 10.087 Combi Clamp Standard
REF. 90.021 for rails 25 x 10 mm
REF. 90.022 for rails 35 x 8 mm
REF. 90.026 for „T“ rails 35 x 8 mm
REF. 90.027 for „T“ rails 35 x 10 mm

for ARGUS Docking Stations
REF 90.041 for rails 25 x 10 mm



Rail Upgrade KIT

additional equipment for the Combi Clamp Standard

REF. 10.108 for rails 25 x 10 mm
REF. 10.109 for rails 35 x 8 mm
REF. 10.110 for „T“ rails 35 x 8 mm
REF. 10.111 for „T“ rails 35 x 10 mm



Interface cable

compatible with ARGUS pumps

- interface RS-232
- updates, configurations, read outs
- software packages on request

REF. 10.093

Cables for docking stations

- interface RS-232 for PC and PDMS connection
- interface RS-242 cascade solution for docking stations



REF. 12.041 Interface RS-232
REF. 12.046 Cascade Solution



Plugs for staff alert systems

compatible with ARGUS pumps

- cable length 2 m

REF. 12.042 for infusion pumps

REF. 10.059 for syringe pumps



DC plug

compatible with ARGUS infusion and syringe pumps

- 20 V
- cable length 3 m

REF. 10.131



Pressure gauge with stopcock

- calibration tool

REF. 12.035



Syringe calibration components

to be used for ARGUS syringe pumps

- part 1: dia. 31 mm, length 20 mm
- part 3: dia. 17 mm, length 120 mm

REF. 11.194 / Part 1

REF. 10.153 / Part 3

We care for every drop



CODAN ARGUS AG, Aarestrasse 13, P.O. Box, CH-3627 Heimberg
Phone +41 33 438 13 38, Fax +41 33 437 11 41, info@codanargus.com
www.codanargus.com

